


# ΜΕΤΡΗΣΗ ΠΡΩΤΕΪΝΩΝ ΜΕ ΤΗΝ ΜΕΘΟΔΟ MULTIPLEX ELISA

Εμβιομηχανική και Βιοϊατρική Τεχνολογία  
Τμήμα Μηχανολόγων Μηχανικών | Ε.Μ.Π.

Χειμερινό Εξάμηνο 2015



# Lab Objective

- Utilize the bead-based multiplex ELISA method to quantify the concentration of Albumin protein in urine samples



# Intro: ELISA

- A standard biochemical method for quantifying the concentration of a protein in solution
  - Utilize 2 antibodies (capture ab, detection ab)



- Takes place in 96-well plates, where capture ab is immobilized
- Quantifies 1 protein per solution
- Utilizes a standard curve to calibrate detected signal


# Albumin

- Main component protein of human blood plasma
- Regulates the pressure of blood
- albuminuria: having too much albumin in the urine, if persistent it is linked to Chronic Kidney Disease (CKD)
  - 30 to 300 mg of albumin means albuminuria

# Bead-based Multiplex ELISA


- Bead based assay

- Immobilize capture ab on beads rather than well walls
- Then add beads to sample to be quantified



- Multiplex assay

- A type of assay that simultaneously measures multiple analytes (proteins) in a single run.
- Here: add to each sample several kinds of beads, each kind can quantify a different protein
- Luminex xMAP technology is one such example



# xMAP Multiplex ELISA Technology

1. Color-coded beads, pre-coated with analyte-specific capture antibody for the molecule of interest are added.
  2. Analyte-specific antibodies capture the analyte of interest.
  3. Biotinylated detection antibodies specific to the analyte of interest are added and form an antibody-antigen sandwich.
  4. Phycoerythrin (PE)-conjugated Streptavidin is added.
  5. The beads are read by a dual-laser flow-based detection instrument which excites the internal dyes marking the beads set and a second laser excites PE, the fluorescent dye on the reporter molecule.
- The system is capable of measuring potentially up to 100 analytes simultaneously in a small sample volume (25–50  $\mu\text{L}$ ).

# Advantages of Multiplex ELISA

- **Speed/High Throughput:** each microsphere serves as an individual test, thus a large number of different assays can be performed and analyzed simultaneously.
- **Accuracy:** generates real-time analysis and accurate quantification of the biological interactions, limited sample handling, and decreased time and cost.

# Standard vs Multiplex Bead-Based ELISA

- Standard ELISA
  - Quantify just one analyte at a time → slow, low throughput
  - Widely utilized all over the world
- MBAA techniques
  - Can be multiplexed → can quantify multiple proteins per sample → high throughput
  - May suffer from antibody cross-reactivity that causes distorted results
    - Only high-quality antibodies can be used
  - Require specialized instruments