

Δημιουργία Οπτικού Σήματος Σε Βιολογικά Δείγματα Μέσω Φθορισμού


Δημήτριος Τζεράνης, Ph.D.

Εμβιομηχανική και Βιοϊατρική Τεχνολογία
Τμήμα Μηχανολόγων Μηχανικών | Ε.Μ.Π.

Χειμερινό Εξάμηνο 2015

The Three Components of Imaging

- Three major components



Emission generation

- Fluorescence

Optical system

- Microscope
- Endoscope


Emission detection

- Detector & sampling

Presentation Overview


- Introduction to fluorescence
- Methods of fluorescent staining
- Förster sensors
- Key Applications of Fluorescence

Intro to Fluorescence



Fluorescence


- An optical phenomenon, where a molecule (fluorophore) absorbs a photon, gets excited, and shortly later emits a photon of less energy and returns to ground state



- Key properties of each fluorophore
 - Absorption spectrum
 - Emission spectrum
 - Emission lifetime


Fluorescence

- Fluorescence absorption and emission is described by a Jablowski energy diagram




Absorption and Emission Spectra

- Absorption spectrum $\alpha(\lambda)$: describes the ability of light of wavelength λ to excite the fluorophore
- Emission spectrum $\varepsilon(\lambda)$: describes the probability that a photon emitted by the fluorophore will have wavelength λ
 - $\varepsilon(\lambda)$ is always red-shifted compared to $\alpha(\lambda)$



Absorption and Emission Spectra

- The spectra $\alpha(\lambda)$, $\varepsilon(\lambda)$ are properties of each fluorophore




Fluorescence Lifetime


- There is a finite time delay between photon absorbance and photon emission



- This delay follows an exponential probability distribution with parameter τ : lifetime
 - Characteristic of a fluorophore
 - On the order of 1 nsec




Fluorescence Staining




Overview: Fluorescence in Biological & Medical Imaging

- Fluorescent molecules (fluorophores) can be visualized using microscopes (see next presentation)



- Utilize fluorophores to label structures or proteins of interests *in vitro* or *in vivo*




Types of Fluorophores Utilized in Biological & Medical Imaging

- Tissues/cells contain molecules that are fluorescent
 - Autofluorescence
- Usually scientists add molecules that fluorescently stain/tag a molecule/structure of interest
 - Fluorescent Stains
 - Fluorescently-conjugated antibodies (immunofluorescence)
 - Fluorescent proteins
 - Fluorescent labeling of peptide tags

Intrinsic Fluorescence of Tissues & Cells (Auto-fluorescence)

- Cells/tissues contain several weakly-fluorescent molecules




Key intrinsic fluorophores found in tissues
(Zipfel et al. 2003)

Imaging skin *in vivo* using intrinsic
fluorophore (Radosevich et al.2008)

Fluorescent Stains

- Chemical compounds that bind specifically to particular molecules found in tissues
 - Nucleic Acid Stains (e.g. DAPI)
 - Actin Stains (phalloidin)
 - Microtubules (paclitaxel)
 - Mitochondria (mitotracker)
 - Cytosol (cell trackers)
 - Cell membrane (WGA)
 - Live VS Dead cells





3T3 cells stained with DAPI (blue), phalloidin (green), and mitotracker (red) (Nikon)

- Cells add the dye to sample → dye diffuses and binds target
- Applied extensively in cell imaging (*in vitro*) or *ex vivo* tissue imaging (thin sections)

Immunofluorescence (IF)


- Utilize a fluorescently-conjugated secondary antibody
 - Binds primary antibody that has bound target protein of interest
- Stain particular proteins in fixed cells *in vitro* or thin tissue sections *ex vivo*
- Can use multiple (usually up to 3) sec. antibodies conjugated to different fluorophores




IF staining of vinculin (red), β -actin (green), and α SMA (blue) in contractile cells. Hinz et al. 2007

Fluorescent Proteins (FP)

- Intrinsically fluorescent proteins, aprox. 25 kDa size
- Green fluorescent protein (GFP) first discovered in *Aequorea victoria* in the 1960s
- Multiple mutants have developed with different properties
 - Emission spectrum (peak), monomeric, stability





GFP domain
(Zhang, 2002)




Fluorescent Labeling via FPs

- Most proteins are not fluorescent
- Can make a protein fluorescent by fusing it with a FP
- Express in live cells *in vitro* or animals


Protein



<http://www.nanobme.org/>




<http://en.dogeno.us/>




Examples of Application of FPs in Animals and Plants

- Generate transgenic animals/plants that express their native proteins fused to a FP
 - Need to incorporate the protein-FP fusion in the animal DNA
 - Enables *in vivo* imaging of particular markers!



Studying neurodegenerative diseases using GFP-mice and *in vivo* imaging (Misgeld et al. 2006)



Imaging fimbrin-GFP fusion in live Arabidopsis roots (<http://www.noble.org/>)

Nobel Prize in Chemistry 2008

The Nobel Prize in Chemistry 2008



Photo: U. Montan
Osamu Shimomura
Prize share: 1/3



Photo: U. Montan
Martin Chalfie
Prize share: 1/3




Photo: U. Montan
Roger Y. Tsien
Prize share: 1/3

The Nobel Prize in Chemistry 2008 was awarded jointly to Osamu Shimomura, Martin Chalfie and Roger Y. Tsien *"for the discovery and development of the green fluorescent protein, GFP"*.




Fluorescent Labelling via Small Peptide Tags

- Fusing a FP with a protein may affect its function
- Alternative tag a protein with a small peptide, that becomes fluorescent upon binding to a small molecule
 - Examples: FIAsh – tetracysteine system




Fluorescent staining of proteins using the FIAsh dye (Madani et al. 2009; Crivat and Taraska 2012)

FRET Sensors




FRET Principle

- Förster resonance energy transfer (FRET) refers to the energy transfer from a “donor” to an “acceptor” fluorophore when they are located <7 nm away
 - Condition for FRET: donor $\epsilon(\lambda)$ overlaps significantly with acceptor $\alpha(\lambda)$





Normally, a blue laser excites CFP, but not YFP



When CFP, YFP are less than 7nm away
→ Blue laser excites YFP through CFP

FRET Sensors


- Sensors refer to artificially design fusion proteins where FRET takes place in ways that provide us information
 - FRET signal depends on protein conformation
 - Protein conformation depends on the entity we want to quantify
- FRET sensors can be expressed in live cells
 - Provide quantitative data *in vivo*




FRET sensor for quantifying the activity of Src small GTPase(Wang et al. 2005)

FRET Sensors


- Sensors designed so far can quantify:
 - Small GTPases
 - Kinase activity
 - Secondary messengers (Ca^{2+})
 - pH, cation concentration



Ras FRET sensor (Yasuda 2006)




Ca⁺⁺ FRET sensor (Ting2001)




Kinase A FRET sensor (Ting2001)

Types of Samples in Biological Imaging


- Cells *in vitro*



- Tissue sections *ex vivo*




- Animal imaging *in vivo*




Whole Animal Imaging

- Utilize fluorophores that emit in the IR
 - Inject/express in live animals
 - Excite using IR light
- Tomographic detection → computation
→ estimate the emission origin inside animal
 - Prof. V. Ntziachristos (Harvard U. , TUM)
- Utilized experimentally to study tumor progression




Major Applications of Fluorescence




Fluorescence: Utilized in Detection

- Fluorescence provides a cheap and versatile way to detect the amount and spatial distribution of molecules that otherwise provide no optical signal...
 - Label molecule of interest with a fluorophore via the abovementioned methods (dyes, antibodies, proteins, tags)
 - Detect fluorophore via fluorescence → detect molecule of interest



Practical Aspects of Fluorescence Staining

- **Photobleaching** phenomenon : repetitive excitation of fluorophores leads to their eventual destruction
 - Prevent by: reducing illumination exposure, prevent O₂ diffusion




Photobleaching example
(micro.magnet.fsu.edu)

- **Photodamage and phototoxicity**: once activated, fluorophores may react with their environment → cause damage to cells
 - Reduce illumination, use necessary staining

Key Instruments That Utilize Fluorescence in Biology & Medicine

- **Fluorescence platereader**

- Cheap, widely utilized
- Quantify solutions or cells




Biotek.de



- **Flow cytometer**

- Quantify fluorophores in cells or beads



- **Fluorescence microscope**

- see next...

Summary

- Fluorophores absorb photons of particular wavelength and within ~ 1 nsec emit photons of less energy
- Tissues/cells contain several weakly fluorescent molecules (autofluorescence)
- Scientists label molecules of interest in cells/tissues using i) fluorescent stains, ii) fluorophore-conjugated antibodies, iii) fluorescent proteins
- Fluorescent proteins fused with proteins of interest in live cells and animals revolutionized science
- The FRET phenomenon has been utilized to design “sensors” that can quantify the activity or concentration of various molecules *in vivo*
- Fluorescence labeling is utilized extensively to detect molecules that otherwise don't emit optical signal